

Orthomation Network Install

Server – Network Server

1. **Install:** Install OrthoMation on Server from Orthomation CD. Choose Network Server as the install type and except all the defaults. (this does not have to be a true server – just a computer you have designated as the computer to store the media for OrthoMation)
2. **Sharing:** Share the Orthomtn folder on the local C drive. (make sure you share with full read/write access)
3. **Net Directory:** Go to Start, Programs, OrthoMation, Support/Tools, and OrthoMation Utilities. Change the net directory to C:\Orthomtn by using the browse button. Click Set and Yes.
4. OrthoMation will now ask for a key which you can obtain from the GAC TechnoCenter (888-422-2376).

Workstations – Network Workstation

(These instructions apply to each workstation OrthoMation is installed on)

1. **Mapping:** Map to the C:\Orthomtn folder you shared on the OrthoMation Server as O:\. (make sure reconnect at log on is checked)
2. **Install:** Install OrthoMation on the workstation from OrthoMation CD. Choose Network Workstation as the install type and then browse to the O:\ folder when prompted for the User Media directory. Finish the install and then reset the BDE.
3. **Resetting the BDE Net Directory:** Go to Start, Programs, OrthoMation, Tools, and OrthoMation Utilities. Change the Net directory to O:\ by using the browse button. Click Set and Yes.
4. **User Folder:** This setting is done in the OrthoMation Utilities window as well. Verify the user folder is O:\. If it is not then use the browse button to select this folder. Click Sync and Yes.
5. **Display File:** This setting is done in the OrthoMation Utilities window as well. Verify the Display File is O:\dispengl.db. If it is not then use the browse button to select this file. Click Set and Yes.
6. OrthoMation will now ask for a key which you can obtain from the GAC TechnoCenter (888-422-2376).

